

CARL Basic

Set Up

1 Unbox your CARL Basic



Your CARL Ears
& Packing Slip



Your CARL Basic



2 "Click-in" Your First Set of Ears



Note: Your CARL Basic comes with two sets of ears - one Realistic (coloured), and one Durable (translucent). Your ear selection can be used to better suit different procedures.

CARL Basic

Considerations for Success

Be Sure to Select the Correct Ears for your Task!

Otoscopy



Hearing Aid Fitting



Impressions



Cerumen Management



Realistic Ears (Coloured)

Durable Ears (Translucent)

How is Best to Set Up Your CARL?

CARL On-the-Go

Pre-installed Suction Cups on Bottom of CARL for Securing to Surfaces




Arm*

Secure CARL using VESA mount to a monitor arm



Stand*

Secure CARL to speaker stand with suction cups



Considerations for your mounting preference:

- **WHO** needs access to the CARL?
- **WHAT** equipment will we use with CARL and needs to be nearby?
- **HOW** OFTEN will CARL need to change locations?

***Available for Purchase**

Scan the QR code on the back of CARL to go right to our extensive CARL resources at aheadsimulations.com!